

Control - X (x-ray)
 VXRue (DR plate)

LABEL

Shielded

NP

Anatomy	Projection	kVp	mAs
Chest	PA		
	Lateral		
Ribs/Chest	AP/PA	80	8
	Oblique		
	Lower		
Sternum	PA Oblique		
	Lateral		
Soft Tissue Neck	AP		
	Lateral		
Abdomen	AP	78	18
Cervical Spine	Lateral	80	8
	Odontoid		
	AP/Oblique	78	6.3
Thoracic Spine	AP	74	12.5
	Lateral	80	20
	Swimmers		
Lumbar Spine	AP/Upright/Obl	78	8
	Lateral	80	10
	L5-S1 Spot		
Sacrum/Coccyx	AP		
	Lateral		
SI Joints	AP/ Oblique		
Pelvis	AP		
Hips	AP/ Lateral		
Femur	Proximal AP/Lateral	76	6.3
	Distal AP/Lateral		
Knee	AP/Lat/Obl		
	Tunnel		
	Skyline		
Tib/Fib	AP/Lateral	68	2
Ankle	AP/Oblique		
	Lateral		

Full view

Anatomy	Projection	kVp	mAs
Foot	AP/Oblique	58	1.6
	Lateral	60	1.6
Toes	AP/Obl/Lat		
Calcaneus	Axial		
	Lateral		
Shoulder	AP		
	Oblique		
	Axial		
AC Joint	AP		
Clavicle	AP/Axial		
SC Joints	AP/PA		
	Oblique		
Scapula	AP		
	Lateral		
Humerus	AP	80	6.3
	Lateral		
Forearm	AP/Lateral		
Elbow	AP/Lat/Obl		
Hand/Wrist	PA/Oblique	60	1.8
	Lateral		
Scaphoid	PA/Obl/Lat		
Finger	PA/Obl/Lat		
Skull	Townes		
	PA AP	88	16
	Lateral	78	12.5
Facial Bones	PA Waters		
	Caldwell		
	Lateral		
Nasal Bones	PA Waters		
	Lateral		
TMJ	Townes		
	Oblique		
Mandible	AP/PA		
	Oblique		
	Lateral		
	Towns		
Orbits	PA Waters		
	Lateral		